



Compact Compressor Nebulizer

Item #	Product Description	UOM
MQ5800	Compact Compressor Nebulizer with Disposable Neb Kit	1/ea

- Dishwasher Safe NebJet.
- Lightweight, quiet, easy-to-Use.
- Durable, Piston-Driven Compressor.
- Integrated Handle for transport.
- Compact, Portable Design.

specifications

Liter flow range	4 to 8L/min.
Nebulizer Medication Capacity	5mL
Particle Size	0.5 µm to 5.0 µm
Compressor Pressure Range	30 to 36 psi
Operating Pressure Range	8 to 16 psi
Operating Temperature Range	10° to 40°C (50° to 104°F)
Operating Humidity Range	10% to 95% RH
Storage Temperature Range	Storage Temperature Range
Storage Humidity Range	10% to 95% RH
Power Requirements	120V, 50-60 Hz
Power Consumption	80 VA
Noise Level	55 dBA
Weight	4.3 lb
Warranty	Limited 5 year

accessories

NEB KIT 500	Disposable Neb Kit with Mouthpiece, "T," and 7' Tubing	50/cs
MASK 001A	Adult Aerosol Mask	50/cs
18254	Reusable Neb Kit with Mouthpiece, "T," and 7' Tubing	10/cs