



HPCPS
E0570
A7005



Panda



Pink Panda

Panda Pediatric Compressor Nebulizer

Item #	Product Description	UOM
MQ6003	Panda Pediatric Nebulizer with Disposable Neb Kit	1/ea
MQ6003R	Panda Pediatric Nebulizer with Reusable & Disposable Kit	1/ea
MQ6004	Panda Pediatric Nebulizer, with Disposable Neb Kit, no carry bag	1/ea
MQ6005	Panda Pediatric Nebulizer, Pink with Disposable Neb Kit	1/ea
MQ6005R	Panda Pediatric Nebulizer, Pink with Reusable & Disposable Neb Kit	1/ea

specifications

Liter Flow Range	4 to 8L/min.
Nebulizer Medication Capacity	5mL
Particle Size	0.5 µm to 5.0 µm
Compressor Pressure Range	30 to 36 psi
Operating Pressure Range	8 to 16 psi
Operating Temperature Range	10° to 40°C (50° to 104°F)
Operating Humidity Range	10% to 95% RH
Storage Temperature Range	-20° to 70°C (-4° to 158° F)
Storage Humidity Range	10% to 95% RH
Power Requirements	120V, 50-60 Hz
Power Consumption	80 VA
Noise Level	55 dBA
Weight	4.3 lb
Warranty	Limited 5 Year

- Child Friendly Design.
- 6003 and 6005 models include disposable neb kit, pediatric mask, tubing, and a carry bag.(6004 models exclude carry bag).
- Durable, piston driven system.

accessories

NEB KIT 500	Disposable Neb Kit with Mouthpiece, "T," and 7' Tubing	50/cs
MASK 001P	Pediatric Aerosol Mask	50/cs
18254	Reusable Neb Kit with Mouthpiece, "T," and 7' Tubing	10/cs