

drive MEDQUIP

NEW  
\*\*\*\*\*1P  
ACKHCPCS  
A4614NEW  
\*\*\*\*\*HCPCS  
A7005

## COMPARE TO

PARI®

RESPIRONICS®

022F81

HS870



## Peak Flow Meter

## Medel® JetNeb Plus

Item #	Product Description	UOM
MQ8500	Peak Flow Meter	1/ea

- Suitable for adults and children.
- Easy-to-read.
- Three zone meter (green, yellow, red).
- Lightweight, portable, pocket-size.
- Convenient Kickstand handle.
- Easy to clean.

## specifications

Dimensions	138 x 45 x 35mm
Weight	42 g
Measurement Range	60-720 L/min (ATS) 60-800 L/min
Accuracy	+/- 10% or 10 L/min
Repeatability	<5 L/min
Warranty	Limited 1 Year

Item #	Product Description	UOM
18260	JetNeb Plus	12/cs

- The Medel® JetNeb Plus breath enhanced reusable nebulizer utilizes a dual valve system improving clinical efficacy and treatment efficiency.
- Through superior product engineering, the Medel® JetNeb Plus uses an inspiratory valve that delivers more medication while the patient breathes in, and an expiratory valve that wastes less when the patient breathes out.
- Optimized average particle size of 3.4 µm for targeted deep lung therapy.
- Up to 92.6% of medication is aerosolized into fine particle mist.

## specifications

Flow	5.1 L/min
Output	0.45 ml
Output Range	0.02 ml / min
Nebulization Rate	0.3 ml / min
Average Particle Size	3.4 µm
Fraction of Fine Particles	92.6%
Maximum Filling Volume	8 ml
Minimum Filling Volume	2 ml