



drive MEDQUIP



HPCPS
E0603



HPCPS
E0603



Dual Channel Breast Pump

Item #	Product Description	UOM
MQ9100	Dual Channel Breast Pump	1/ea

- LCD Display provides status of programs and functions.
- Adjustable Multi-level Dual Channel for Optimum Expression.
- Includes Two Soft Silicone Cushion Breast Cups.
- Anti-backflow System Prevents Milk from Entering Motor.
- Bottle Feeding Adapter Kit & Storage Containers with Lids.
- Alarm Function Reminder for next Pumping Session.
- Bisphenol A (BPA) Free.
- May be powered with four AA Alkaline Batteries.
- A/C Power Supply Included.
- Convenient Carrying Bag.

specifications

Warranty	One Year Limited Consumer Warranty
----------	------------------------------------

Dual Channel Breast Pump

Item #	Product Description	UOM
MQ9120	Dual Channel Breast Pump	1/ea

- LED Display provides status of programs and functions.
- Adjustable Multi-level Dual Channel for Optimum Expression.
- Includes Two Soft Silicone Cushion Breast Cups.
- Bottle Feeding Adapter Kit & Storage Containers with Lids.
- Closed System Anti-backflow Diaphragm prevents milk from entering the tubing or electric motor.
- Bisphenol A (BPA) Free.
- May be powered with four AA Alkaline Batteries.
- A/C Power Supply Included.
- Convenient Carrying Tote.
- One Year Limited Consumer Warranty.

specifications

Power	4 AA Battery/AC
Type	Double
Display	LED
Vacuum Chamber	2
Suction Levels	7
Vacuum Strength	17mmHg ~ 250mmHg
Speed Setting	4 Preset Speeds
Standard Accessories	2 Bottle Caps 2 Silicone Cushions