

MSS Tilt & Recline

Available
1 Pack

Model #:
Sold 1 per box

MS 3000N

MSS Tilt & Recline
Seating System



The MSS Tilt & Recline seating system is designed for children with moderate to severe physical involvement. The newly designed chair has a hydraulic actuator that activates the tilt in space feature. One model now provides the same growth range as our previous two models. The basic model includes all the features below.

Weight Capacity: 100 lbs

Warranty: Limited Lifetime

Features

- Multi-Axis Headrest—Height, depth and angle adjustable. It can be placed off center.
- Seat—Depth Adjustable.
- Back—Height and angle adjustable.
- Harness—Adjustable zippered "H" harness for anterior chest support.
- Pelvic Belt—Adjustable in length.
- Laterals—Height and width adjustable. They assist in stabilizing and aligning trunk and pelvis. Each lateral is individually adjustable to control unwanted weight shifts.
- * Armrest—Height and angle adjustable.
- Adductors—Two sizes are included to facilitate a range of chair widths. The adductors control lower extremity positions and prevent adduction of user's legs.
- Abductor—Height and depth adjustable.
- Footrests—Height, depth and angle adjustable. They include foot sandals with toe loop. Each footrest is individually adjustable to custom fit.
- * Removable and washable.
- Limited Lifetime Warranty.

Specifications

Item #	Seat Width	Seat Depth	Back Height	Overall Width	Seat to Footrest Height	Seat to Floor Height	Distance Between Laterals	Headrest Height Adds
MS 3000N	9"-10.5"	3.5"-12"	14"-19"	18"	4"- 10.5"	14"	9"-15"	4"-10"
Item #	Product Weight		Carton Shipping		Weight Limit Weight			
MS 3000N	35 lbs.		40 lbs.		100 lbs.			

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.