

# MSS Tilt & Recline Accessories



MS 0014N

## Base

**Model #:**  
Sold 1 per box



**MS 001 4N**  
**MS 001 5N**

Hi Base for MS 3000N  
Low Base for MS 3000N

Wheeled mobility bases allow the MSS Tilt & Recline to be used as a floor sitter or high chair that is brought up to table top height.

**Weight Capacity:**

100 lbs

**Warranty:**

Limited Lifetime



MS 0015N

## Specifications

Item #	Seat to Floor Height (chair mounted on base)	Product Weight	Carton Shipping Weight	Weight Limit
MS 001 4N	24"	6 lbs.	9 lbs.	100 lbs.
MS 001 5N	17"	6 lbs.	9 lbs.	100 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.

**IP = 1-Pack Available**

To order call toll free: 877.224.0946

WENZELITE

Wenzel<sup>®</sup> Lite Recline  
**drive**