

# MSS Tilt & Recline

The MSS Tilt & Recline is designed for children with moderate to severe physical involvement. It features an easy tilt and recline mechanism with up to 45 degrees tilt and the back reclines from 80° up to 140°. The basic model includes all the features listed below.

MS 2000 Small, 1/cs

MS 3000 Medium, 1/cs

## FEATURES:

- Footrest – adjustable in height and depth.
- Seat – adjustable in depth.
- Armrest – adjustable in height and angle.
- Back – adjustable in height and angle.
- Headrest – adjustable in height, depth and angle. It rotates and can be placed off center.
- Lateral supports – individually adjustable in height and width.
- Abductor – adjustable in height and depth. It can be placed between the thighs or turned around and placed between the knees.
- Belt – adjustable pelvic belt.
- Tray – horizontal positioning device.
- Harness – zippered trunk harness in an H configuration for anterior chest support.
- Cushion covers are removable and washable.
- Two year limited warranty.

## BASE:

MS 0014 Hi Base for MSS Tilt & Recline, 1/cs

MS 0015 Low Base for MSS Tilt & Recline, 1/cs

Wheeled mobility bases fit both sizes MSS Tilt & Recline chairs. It allows the MSS Tilt & Recline to be used as floor sitter or a high chair. The Hi Base will bring the seating system up to table height.

DIMENSIONS	MS 2000	MS3000
SEAT WIDTH	9"	10.5"
SEAT DEPTH	5" - 14"	5" - 14"
BACK HEIGHT	10" - 12.5"	12.5" - 16"
SEAT TO FOOTREST HEIGHT	5" - 10.5"	5" - 10.5"
FOOTREST DEPTH	7.5"	7.5"
ARMREST TO SEAT HEIGHT	5" - 9"	5" - 9"
FLOOR TO SEAT HEIGHT		
AS FLOOR SITTER	13.5"	13.5"
WITH LOW BASE	16"	16"
WITH HIGH BASE	22.5"	22.5"
SPACE BETWEEN LATERALS	7" - 11"	7" - 11"
HEADREST HEIGHT ADJUSTMENT	Adds 3" - 7"	Adds 3" - 7"
WEIGHT LIMIT	75 lbs.	100 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.



positioning  
seating

Wenzel<sup>®</sup>iteRe/hab

To order call toll free: 877.224.0946

drive  
WWW.DRIVEMEDICAL.COM