



NEBKIT 500



NEBKIT 500R-12



NEBKIT 600

Nebulizer Kit (Disposable)

Nebulizer with T-adapter, 7' tubing, mouthpiece and reservoir tube. Provides 80% of particles in respirable range.

NEB KIT 500
NEB KIT 600

50/cs
Nebulizer Cup,
7' Crush Resistant
Tubing, Adult Aerosol
Mask, 50/cs

NEB KIT 700

Nebulizer Cup,
7' Crush Resistant
Tubing, Pediatric Aerosol
Mask, 50/cs

NEB KIT 500R-12

Retail, 12/bx



NEBKIT 700

COMPARE TO:

	drive	Hudson®	Salter®
NEBKIT 500	1883	8900	
NEBKIT 600	1885	8924	
NEBKIT 700	1886	8906	



COMPARE TO:

	drive	Hudson®	Salter®
18254	1901	8660	



Reusable Nebulizer Kit

18254

Contains Mouthpiece,
"T" reusable Neb Cup
and 7' Tubing, 10/cs

FEATURES:

- Anti-spill "T" design.
- Leak-proof and secure seal with simple 1/4 turn assembly design.
- Very low residual volume (0.7cc).
- High medication bowl capacity (cc)
- Graduated volume (cc) markings on bowl.
- Nebulizes in upright positions or at angles up to 45°.
- Latex free.
- Dishwasher safe.

Humidifier Bottle

HUM 001

For Use with all Oxygen Sources:
Concentrators, Liquid, Tank,
and Wall Outlets, 50/cs

FEATURES:

- Audible alarm activates at 6 psi.
- Rigid and durable wing-nut and lid.
- No aluminum fitting in stem: it eliminates shredded aluminum diffuser so often found in cheaper models. Soft plastic diffuser reduces clogging due to mineral build-up.
- Easy-on/easy-off lid and jar thread, provides positive seal.
- Translucent jar provides improved visibility of jar contents.



HUM 001

COMPARE TO:

	drive	Hudson®	Salter®
HUM 001	N/A	7100	

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.