







Otter Bathing System

		Product	Product	Product	Product				
	Product	Seat	Back	Outside	User	Product	Product Weight	Product Carton	
Item #	Size	Depth	Height	Overall Length	Height	Weight	Capacity	Shipping Weight	UOM
OT 1000	Small	13"	18"	36.5"	Up to 36"	11 lbs.	60 lbs.	14 lbs.	1/bx
OT 2000	Medium	13"	25"	42"	32"- 50"	12 lbs.	120 lbs.	15 lbs.	1/bx
OT 3000	Large	16.5"	32"	50"	46"- 68"	14 lbs.	160 lbs.	18 lbs.	1/bx
OT 1000 SF	Small	13"	18"	36.5"	Up to 36"	11 lbs.	60 lbs.	14 lbs.	1/bx
OT 2000 SF	Medium	13"	25"	42"	32"- 50"	12 lbs.	120 lbs.	15 lbs.	1/bx
OT 3000 SF	Large	16.5"	32"	50"	46"- 68"	14 lbs.	160 lbs.	18 lbs.	1/bx

- Plastic frame with a choice of standard (Figure B) or soft (Figure C) removable and washable fabric.
- Standard Fabric: Vinyl covered nylon dries easily. (Figure B)
- Soft Fabric: Polyester knit fabric includes a trunk strap to assist in lateral support. Features additional padding between the fabric and Otter frame. (Figure C)
- Angle-adjustable seat and back have five angle adjustments at 0°, 22.5°, 45°, 67.5° and 90°.

specifications	
Seat Dimensions	14.5" (W) x 2"- 7" (H)
Overall Outside Width	17"
Warranty	Limited Lifetime

- Seat and back Uni-Bars provide one-handed angle adjustments. (Figure A)
- Adjustments can be made with user seated in chair.
- Height- and width-adjustable head cushions offer support and comfort. They can also be used as lateral supports.
- Removable positioning straps can be placed anywhere along the length of the frame.
- Leg straps control abduction and adduction.
- Adjustable slip-resistant legs raise chair up to 7".
- Folds flat for easy storage.