

PrimeCare Model #P902

LONG TERM CARE BED




Low Bed
Height 7"



76", 80", 84"
Mattress Length



Optional Width
Expansion to 42"



Floor Lock
System



Wall Bumper &
Steer Control



Shown with Optional Bed Ends,
Assist Bar and Embedded Staff Control

Features

- Easy-to-use Hand Pendant controls positioning for maximum resident comfort and versatility for the caregiver.
- Comfort Chair and Cardiac Chair (with Trendelenburg/Reverse Trendelenburg option) adjustability helps provide added comfort and relaxation positioning for the resident and facilitates caregiver accessibility.
- Edema foot ratchet elevates the foot end of the bed
- Trendelenburg/Reverse Trendelenburg capable positioning with Staff Control Option is ideal for rehab and sub acute care, provides therapeutic benefits that can aid in circulation and ease discomfort.
- Roll at any height convenience, plus steer control provides easy maneuverability
- Easy access Bar-Pedal Floor Lock system at the foot end. Press Bar down to engage floor locks; disengage using release pedals on either side of the caster housing.
- Adjustable mattress retainer and wall saver are standard
- All Steel frame and no sag slat deck construction offers superior strength, durability, mattress support and is easy to clean
- Options: 3 Bed End Styles/5 attractive wood grain Laminate colors; Assist Bars/Rails; Embedded Staff Control w/lock out; Trapeze; Battery backup

Specifications


Sleep Surface	76", 80", 84"
Base Width	35"
Optional Width Expansion	42"
Working Height	30.5"
Super Low Height	7"
Trendelenburg/Reverse Trendelenburg Capable	With Staff Control Option
Weight Capacity	500 lbs.*
Warranty	15 Year on Frame, 3 Year on Electric and Lifetime on All Welds

Head and Footboards

LONG TERM CARE BED FURNITURE



35" HEAD AND FOOT BOARDS

BROWNLEE WOODS RECTANGULAR (BW)





42" HEAD AND FOOT BOARDS



MILL CREEK ARCHED (MC)



SMOKEY HOLLOW (SM)



38" HEAD AND FOOT BOARDS (FITS BARIATRIC BEDS ONLY)



35" HEAD/FOOT BOARDS (PAIR) FITS P301, P601, P902

Finish	Brown Lee Woods	Mill Creek	Smokey Hollow
35" Solar Oak	35BWHF-SO	35MCHF-SO	35SHHF-SO
35" Wild Cherry	35BWHF-WC	35MCHF-WC	35SHHF-WC
35" Montana Walnut	35BWHF-MW	35MCHF-MW	35SHHF-MW
35" Figured Mahogany	35BWHF-FM	35MCHF-FM	35SHHF-FM
35" Windsor Mahogany	35BWHF-WM	35MCHF-WM	35SHHF-WM

42" HEAD/FOOT BOARDS (PAIR) FITS P301, P601, P902

Finish	Brown Lee Woods	Mill Creek	Smokey Hollow
42" Solar Oak	42BWHF-SO	42MCHF-SO	42SHHF-SO
42" Wild Cherry	42BWHF-WC	42MCHF-WC	42SHHF-WC
42" Montana Walnut	42BWHF-MW	42MCHF-MW	42SHHF-MW
42" Figured Mahogany	42BWHF-FM	42MCHF-FM	42SHHF-FM
42" Windsor Mahogany	42BWHF-WM	42MCHF-WM	42SHHF-WM

38" HEAD/FOOT BOARDS (PAIR) FITS P1752, P2002

Finish	38" Sculpted
38" Solar Oak	PBAHF-SO
38" Wild Cherry	PBAHF-WC
38" Montana Walnut	PBAHF-MW
38" Figured Mahogany	PBAHF-FM
38" Windsor Mahogany	PBAHF-WM


See www.drivemedical.com for complete listing of bed ends with staff control options.

AVAILABLE LAMINATE FINISHES


SOLAR OAK




WILD CHERRY




MONTANA WALNUT



FIGURED MAHOGANY



WINDSOR MAHOGANY



The color samples shown are only printed representations of our laminate offering, actual colors may vary.
Laminate samples, Item #WP-SAMP2, are available upon request.

Product configurations and specifications subject to change without notice.



99 Seaview Blvd., Port Washington, NY 11050



800.371. 2266



www.drivemedical.com

Accessories

LONG TERM CARE BED FURNITURE

ASSIST RAILS



P300-07AE



P300-08AE



P902ABE



P902ARE2



PBASR

ITEM #	FITS WITH
P300-07AE	P301, P601 Bed Frame
P300-08AE	P301, P601 Bed Frame
P902ABE	P902 Bed Frame
P902ARE2	P902 Bed Frame
PBASR	P1752, P2002 Bed Frame

WIDTH EXTENSION



P300-04E



P902WE

ITEM #	FITS WITH
P300-04E	P301, P601 Bed Frame
P902WE	P902 Bed Frame