





PC 1000 R



PC 1000 BL

HPCPS
E1399

Pinniped Pediatric Commode

Item #	Color	UOM
PC 1000 R	Red 	1/bx
PC 1000 BL	Blue 	1/bx

- All-in-one welded steel frame available in red and blue.
- Plastic seat and contoured back.
- Durable snap-on toilet seat.
- Height adjustable seat allows child to firmly place their feet on the floor for comfort, safety & security.
- Adjustable hook & loop chest strap for additional support.
- Molded plastic armrests provide additional comfort.
- Rubber non-skid tips ensure that the legs will not mark the floor.
- Complete with 7.5 qt removable commode bucket, metal handle and cover.

specifications

Seat Depth	15"
Seat Width	13"
Seat Height	10.5" – 14.5"
Armrest Height (from floor)	15.5" - 19.5"
Armrest Pad	2" (W) x 9" (D)
Overall Width	17"
Toilet Opening	6.5" W x 8" D
Weight	9.5 lbs
Carton Shipping Weight	13 lbs
Weight Limit	250 lbs
Warranty	Limited Lifetime

accessories

RTL12085	Biodegradable Sanitary Commode Liner	12 liners/box, 6bx/cs
----------	--------------------------------------	-----------------------