

Pediatric Anterior Safety Roller

PE TYKE	Tyke, Red, 1/cs
PE 1000 BL	Pediatric, Blue, 1/cs
PE 1000 WB	Pediatric (Extra Wide), Steel, Burgundy, 1/cs
PE 1000 WBK	Pediatric (Extra Wide), Steel, Black, 1/cs
PE 1000 XL	Pediatric (Adj. Width), Steel, 1/cs

The Pediatric Anterior Safety Roller is available in four sizes to accommodate children from early childhood through young adulthood. The Safety Roller's adjustments and compatible modular components provide substantial built-in growth. This in conjunction with its sturdy and functional design allow for many years of use. The Safety Roller builds confidence in children as they encourage safe ambulation and a sense of freedom. The child is never without support and is always in full control.

FEATURES:

- Applying pressure to the handlebar activates braking mechanism. When brake is engaged, the brake stopper hits the wheel and the wheels cannot roll. Ideal for patients with limited hand function and/or limited cognizance. The Safety Roller halts in step with the patient and the braking mechanism automatically engages when patient stumbles.
- Height adjustable in 1" increments.
- Adjustable brake spring tension.
- Patient is situated in the center of the Safety Roller, which provides maximum support and improved posture.
- Skid resistant, non-marring, easy rolling neoprene wheels. Ideal for indoor and outdoor use.
- Folds easily and stands on its own in the folded position.
- Limited Lifetime Warranty.

NEW!
NOW AVAILABLE IN
LIGHTWEIGHT ALUMINUM!



PE 1000 BL



PE 1000 XL



PE 1000 BL

ITEM #	HANDLEBAR HEIGHT RANGE	OVERALL WIDTH	WIDTH INSIDE HANDLE GRIP	DEPTH OPENED AT BASE	WEIGHT LIMIT	PRODUCT WEIGHT	CARTON SHIPPING WEIGHT
PE TYKE	15.5"-18.5"	19"	14"	25"	100 lbs.	14 lbs.	19 lbs.
PE 1000 BL	22" - 28"	20"	15"	25"	200 lbs.	10 lbs.	15 lbs.
PE 1000 WB/BK	24" - 30"	21.25"	17"	25"	200 lbs.	16 lbs.	21 lbs.
PE 1000 XL	23" - 29"	22.5"-27.5"	17.5"-22.5"	25"	400 lbs.	18 lbs.	23 lbs.

* For a complete line of accessories, see pages 426-430.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.