



Drive Medical is proud to support the American Cancer Society's efforts to finish the fight against breast cancer through awareness, early detection, and research.



MAKING STRIDES
Against Breast Cancer®



Aluminum Rollator, 6" Casters

Fold-up and Removable Back Support, Universal Seat Height, Padded Seat, Loop Locks

Item #	Color	UOM
RTL10261BC	Pink 	1/cs

- Universal height adjustment on frame allows seat height to be adjusted from 18" to 22" in 1" increments.
- Comes with seamless padded seat with zippered pouch under seat for added privacy and security of personal items.
- Removable, hinged padded backrest can be folded up or down.
- Caster journals are offset to provide strength.
- Easy to use deluxe loop locks.
- Brakes with serrated edges provide firm hold.
- Ergonomic handles are easy to grip and are adjustable in height with self-threading knob.

specifications	
Handle Height	29.5" - 38"
Unit Dimensions	22" (W) x 27" (L)
Seat Dimensions	12" (D) x 13" (W) x 18" - 22" (H)
Weight	14 lbs.
Carton Shipping Weight	16 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

For more information, call Customer Service at 800-371-2266, contact your Sales Representative, or visit www.drivemedical.com TODAY!