

SIDE POUCHES
MAKE IT EASY
TO CARRY
ITEMS
WITH
YOU!

COMES
WITH
GLIDE
SKIS!

COMES
WITH
5" WHEELS!

EASY FOLDING
DESIGN USES
ELASTIC CORDS!



Deluxe Folding Walker

10263KDR

Deluxe, Folding Travel Walker
with 5" Wheels and Fold Up Legs, 1/cs

FEATURES:


- Easily folds tool free to 1/2 the size of a traditional walker and can be stored in a "Ready Set Go" pouch provided with each walker.
- Easily unfolded and height adjusts tool free.
- Comes with 5" front wheels and rear glide skis that are ideal for most surfaces.
- Two side storage pouches with large pockets for personal items are held in place by Velcro tabs to both side frames.
- Sturdy 1" diameter anodized, extruded aluminum construction ensures maximum strength while remaining lightweight.
- Limited Lifetime Warranty.



LEGS
CONVENIENTLY
FOLD UP
FOR EASY
STORAGE!




WALKER
CAN BE
FOLDED FLAT
FOR STORAGE
& TRANSPORT!



LEGS ARE
ATTACHED WITH
STRONG ELASTIC
CORDS SO LEGS
CAN BE FOLDED!



COMES WITH
READY SET
GO BAG!



PIN MUST GO INTO LOCK COLLAR
BEFORE USE.

SPECIFICATIONS

Item #	Width (inside hand grip)	Width (inside back legs)	Overall Width	Depth (opened at base)	Height (highest setting)	Height (lowest setting)	Product Weight	Carton Shipping Weight	Weight Limit
10263KDR	17"	22"	26.5"	20.5"	37"	29"	7.9 lbs.	9.7 lbs.	300 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

walkers

drive

To order call toll free: 877.224.0946

WWW.DRIVEMEDICAL.COM