



RTL10266-H



RTL10266BK-H



RTL10266WT-H




Hemi Height Ideal For Users 5'2" and Below




A




B



C



D



E




F



Nitro Aluminum Rollator, Hemi Height, 10" Casters

4 Wheel Rollator

Item #	Color	UOM
RTL10266-H	Red 	1/ea
RTL10266BK-H	Black 	1/ea
RTL10266WT-H	White 	1/ea

- Hemi-height ideal for users 5'2" and below.
- Attractive, Euro-style design.
- Brake cable inside frame for added safety.
- Handle height easily adjusts with unique push button. (Figure A)
- Back support height easily adjusts with tool-free thumb screw. (Figure B)
- Caster fork design enhances turning radius.
- Large 10" front casters allow optimal steering and rolling comfort.
- Easily folds with one hand to ultra-compact size for storage. (Figure C)
- Plastic clip holds the Nitro closed when frame is folded. (Figure D)
- Cross-brace design allows for side-to-side folding and added stability.
- Seat is durable and comfortable. (Figure E)
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded. (Figure F)
- Lightweight aluminum frame.

specifications

Handle Height	28" - 31"
Unit Dimensions	23" (W) x 27.75" (L)
Seat Dimensions	18" (W) x 18" (H) x 8.5" (D)
Seat Height	18"
Weight	17 lbs.
Carton Shipping Weight	29 lbs.
Weight Limit	300 lbs.
Warranty (Frame)	Limited Lifetime
Warranty (Brake and Handle Housing)	5 Year Limited