



RTL10266-T



RTL10266BK-T



RTL10266WT-T

Tall Height Ideal For
Users 6'0" and Above.



Nitro Aluminum Rollator, Tall Height, 10" Casters

4 Wheel Rollator

Item #	Color	UOM
RTL10266-T	Red	1/ea
RTL10266BK-T	Black	1/ea
RTL10266WT-T	White	1/ea

- Tall height ideal for users 6'0" and above.
- Attractive, Euro-style design.
- Brake cable inside frame for added safety.
- Handle height easily adjusts with unique push button. (Figure - A)
- Back support height easily adjusts with tool-free thumb screw. (Figure - B)
- Caster fork design enhances turning radius.

- Large 10" front casters allow optimal steering and rolling comfort.
- Easily folds with one hand to ultra-compact size for storage. (Figure - C)
- Plastic clip holds the Nitro closed when frame is folded. (Figure - D)
- Cross-brace design allows for side-to-side folding and added stability.
- Seat is durable and comfortable. (Figure - E)
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when Nitro is open or folded. (Figure - F)
- Lightweight, aluminum frame.

specifications

Handle Height	36" - 41"
Unit Dimensions	23" (W) x 27.75" (L)
Seat Dimensions	18" (W) x 20.5" (H) x 10" (D)
Seat Height	23.6"
Weight	17.5 lbs.
Carton Shipping Weight	30 lbs.
Weight Limit	300 lbs.
Warranty (Frame)	Limited Lifetime
Warranty (Brake and Handle Housing)	5 Year Limited