

Bariatric Aluminum Folding Cane, Height Adjustable

Item #	Color	UOM
RTL10304HD	Black	1/ea

- 500 lb. Weight Capacity.
- Cane folds into four convenient parts for easy storage.
- Attractive wood handle with brass collar.
- Manufactured with sturdy, extruded aluminum tubing.
- Comes with plastic clip to hang cane when folded.
- Handle height adjusts from 33" to 37".

specifications	
Handle Height	33" - 37"
Weight Limit	500 lbs.
Warranty	Limited Lifetime