



























Offset Handle Aluminum Canes Foam Rubber Grip

Item #	Color	UOM
10303-6	Silver	6/cs
10303-1	Silver	1/cs
RTL10306	Black	6/cs
RTL10307	Bronze	6/cs
10308C-6	Silver, Black & Bronze	6/cs

- Ergonomically designed handle with soft Foam Grip provides comfort and security.
- Comes standard with foam rubber grip and wrist strap. (Figure A)
- Manufactured with sturdy, extruded aluminum tubing.
- Easy-to-use, one-button height adjustment with locking ring. Locking ring prevents rattling. (Figure B)
- Available in Silver (main image), Black (Figure C) or Bronze. (Figure D)
- Handle height adjusts from 30" to 39".
- 300 lb. Weight Capacity.

specifications	
Handle Height	30" - 39"
Weight Limit	300 lbs.
Warranty	Limited Lifetime

Offset Handle Cane Combo Pack

Item #	Color	UOM
10309DCW-6	Black Wave, Blue Wave & Bronze Wave	6/cs

- Ergonomically designed handle with soft Foam Grip provides comfort and security.
- Comes standard with foam rubber grip and wrist strap.
- Manufactured with sturdy, extruded aluminum tubing.
- Easy-to-use, one-button height adjustment with locking ring. Locking ring prevents rattling. (Figure A)
- Comes with Black wave (main image), Blue wave (Figure B) and Bronze wave. (Figure C)
- Handle height adjusts from 30" to 39".
- 300 lb. Weight Capacity.

specifications	
Handle Height	30" - 39"
Weight Limit	300 lbs.
Warranty	Limited Lifetime