





RTL10307AT









RTL10306AT







All Terrain Canes

Item #	Color	UOM	
RTL10306AT	Black 🛑	6/cs	
RTL10307AT	Bronze =	6/cs	

- Walk safely and confidently on virtually any surface and under most weather conditions. (Figure A B)
- The retractable tip has four stainless steel spikes to provide a sturdy grip on slippery surfaces.
- No installation required and provides the best traction and stability.
- Other manufacturers' items require separate installation and are less safe.
- To prevent damage to floors, gripping claw should only be used outdoors.
- Manufactured with sturdy, extruded aluminum tubing.
- Available in black or bronze (Figure C D)
- Height adjustable from 30" to 39".
- Retail packaged.

specifications	
Height	30" - 39"
Weight Limit	250 lbs.
Warranty Limited	Lifetime