

FOLDING CANES WITH GEL GRIP HANDLE



**THE NEW GEL GRIP
REDUCES HAND
STRESS AND
FATIGUE!**



Aluminum Folding Cane with Gel Grip, Height Adjustable

10370BK-1	Black, 1/cs
10370BK-6	Black, 6/cs
10370BC-1	Blue Crackle, 1/cs
10370BC-6	Blue Crackle, 6/cs
10370RC-1	Red Crackle, 1/cs
10370RC-6	Red Crackle, 6/cs



**THIS ITEM
IS AVAILABLE
IN A 1-PACK!**

FEATURES:

- Cane folds into 4 convenient parts for easy storage.
- Handle height adjusts in 1" increments from 33" to 37".
- New gel grip reduces hand stress and fatigue.
- **Attractive wood handle .**
- Comes with plastic holster carry case.
- 300 lb. Weight Capacity.
- Limited Lifetime Warranty.



10370RC-1
Red Crackle



10370BC-1
Blue Crackle



10370BK-1
Black

**CANE
DISPLAY
AVAILABLE**

SEE PAGE 102

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

canes & crutches

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM