

K. D. Universal Aluminum Crutches with Accessories

Available **1 Pack** RTL10433 For Individuals 4' 6" to 6'6", 8pr/cs

NEW!

**UNIVERSAL
HEIGHT
ADJUSTMENT
FOR USE BY
CHILD,
ADULT OR
TALL ADULT!**

FEATURES:

- Universal height adjustment for use by a child, adult and tall adult. (individuals 4'6" to 6'6")
- Crutch disassembles, fitting into a retail box 8" x 4.5" x 21.5".
- New Euro – style clips and push pin adjustment make it easy to adjust hand grip.
- Lighter than wood, the aluminum crutch provides stability and durability.
- Non- skid, jumbo-size vinyl contoured tips provide excellent traction.
- Comfortable, durable underarm pad and hand grip.
- Weight capacity: 300 lbs.
- Limited Lifetime Warranty.
- Comes in retail box.



**ATTRACTIVE
RETAIL
PACKAGING!**

**CRUTCH
DISASSEMBLES
FOR
CONVENIENT
STORAGE!**

RTL10433

SPECIFICATIONS

Item #	Height (highest setting) (armrest to floor)	Height (lowest setting) (armrest to floor)	Patient Height (highest)	Patient Height (lowest)	Weight Limit
RTL10433	62"	39"	6'6"	4'6"	300 lbs.

canes & crutches

To order call toll free: 877.224.0946

drive