



I-Walker Aluminum Rollator, 7" Casters

Item #	Color	UOM
RTL10555BL	Blue	1/ea
RTL10555RD	Red	1/ea
RTL10555SL	Silver	1/ea

- Attractive, Euro-style design.
- Handle height easily adjusts with unique push button.
- Back support height easily adjust with tool-free thumb screw.
- Caster fork design enhances turning radius.
- 7" front casters allow excellent steering and rolling comfort.
- Easily folds with one hand to ultra-compact size for storage.
- Plastic clip holds the I-Walker closed when frame is folded.
- Cross-brace design allows for side-to-side folding and added stability.
- Seat is durable and comfortable.
- Handsome and removable zippered storage bag with unique attachment to keep bag securely in place when I-Walker is open or folded.
- Lightweight frame.

specifications	
Handle Height	33.75" - 38.75"
Unit Dimensions	27" (L) x 25.625" (W)
Seat Dimensions	18" (W) x 9.75" (D)
Weight	17 lbs.
Seat Height	20"
Weight Capacity	300 lbs.
Warranty (frame)	Limited Lifetime
Warranty (Brake, Handle & Housing)	5 Years