



11148KDR-1

**PATENT PENDING**



**BOX IS 30% SMALLER THAN CONVENTIONAL COMMODE BOX!**

**350 LB**  
WEIGHT CAPACITY

**Folding Steel Commode**

11148KDR-1 (Retail box, standard seat depth), (Has replaced 11148-1), 1/cs

**FEATURES:**

- Easily opens and folds.
- Folds flat for convenient storage and transportation.
- Durable plastic snap on seat and lid.
- Blue powder coated steel welded construction increases strength and durability.
- Plastic arms for added comfort.
- Comes complete with 12 qt. commode bucket with carry handle, cover and splash shield.
- **Limited Lifetime Warranty.**


**LEGS FOLD UP FOR SMALLEST FOLDED COMMODE AVAILABLE!**



ELASTIC CORD ALLOWS LEG TO BE FOLDED UP TO REDUCE SIZE FOR STORAGE AND TRANSPORT!



**FOLDS FLAT FOR STORAGE!**



EASILY ASSEMBLES WITHOUT TOOLS!

**SPECIFICATIONS**

Item #	Width (between arms)	Outside Width	Width Seat	Depth Seat	Height (highest)	Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
11148KDR-1	18"	22"	13.5"	14.5"	22"	16"	14.8 lbs.	17.8 lbs.	350 lbs.

\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.