



Portable Bath Step

RTL12068

1/cs

FEATURES:

- Makes getting in and out of the bath tub easy and safe.
- Unique modular design allows the multiple steps to be stacked.
- Step height is 4".
- Each additional step adds 2".
- Bath step comes with four removable covers on the surface that locks into another step when stacking.
- Attractive, textured surface is slip- resistant and comfortable.
- 450 lb. weight capacity.
- Dimensions: 19.5" (W) x 16" (D) x 4" (H).
- One Year Limited Warranty.

STEPS CAN BE STACKED, ADDING 2" IN HEIGHT FOR EACH ADDITIONAL STEP.



TOOL FREE

NSTALLATION

PLUGS CAN BE REMOVED WHEN STACKING STEPS TO ENSURE A SECURE FIT.



Plastic Tub Rail

12039

12" White, (Will become RTL12039 in 2011), 1/bx

12044

R1L12039 in 2011), 1/bx 16¹/2" White, (Will become RTL12044 in 2011), 1/bx

FEATURES:

- Attractive, ergonomically designed tub rail provides a comfortable, safe grip.
- Tool free, easy-to-install.
- Width adjusts from 2 3/4" 6 1/2".
- Limited Lifetime Warranty.



SPECIFICATIONS

12044

Item #	Adjustable Clamp (smallest)	Adjustable Clamp (widest)	Product Height	Product Weight	Carton Shipping Weight	Weight Limit
12039	2.75"	6.5"	12"	3.5 lbs.	5.5 lbs.	300 lbs.
12044	2.75"	6.5 "	16.5 "	4 lbs.	6 lbs.	300 lbs.