



## READY SET GO Deluxe Aluminum Shower Chair With Carry Bag

## Shower Chair with Back and Removable Padded Arms

Item #	Product Description	Assembly	UOM
RTL12105KDR	Comes with Carry Bag	Knock Down	2/cs

- Easy, tool-free assembly of legs and back. (Figure A - Figure B)
- Drainage holes in seat and back reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- New pinch-free cover allows for push pins to be depressed without pinching finger. The audible "click" will be heard when secure.
- Legs are height adjustable in 1" increments, and angled with suction style tips to provide additional stability.
- Support collar prevents leg movement and reduces rattle.
- Comes with a durable, nylon carry bag for easy and convenient transport. (Figure C)

specifications	
Outside Legs	19" (W) x 18" (D)
Seat Dimensions	19.75" (W) x 12" (D)
Seat Height	14" - 18"
Weight	8 lbs.
Carton Shipping Weight	19 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

Item #	UOM
12445-2	2/cs
12445-1	1/cs

- Removable padded armrests provide comfort and added security.
- Blow-molded bench and back provide comfort and strength.
- Drainage holes in bench reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- Leg height adjusts in 1" increments.

specifications	
Outside Legs	17.75 (W) x 19.25 (D)
Seat Width	16"
Seat Depth	13", 14.5", 16.25"
Seat Height	16" - 21"
Weight (-1   -2)	9.75 lbs.   19 lbs.
Carton Shipping Weight (-1   -2)	10 lbs.   21 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime