

INTRODUCING OUR.....

Premium Series Bath Bench with Back and Arms



RETAIL PACKAGED!

NEW!



RTL12505 With Back, 4/cs

FEATURES:

- Easy tool free height adjustment.
- Enhanced comfort seat fits all shapes and sizes.
- Tool free removable back, leg and arms.
- Lightweight design has large stable legs with soft non-slip feet.
- 350 lb. weight capacity.



EASY TO USE HEIGHT ADJUSTMENT!

SPECIFICATIONS

Item #	Width (outside legs)	Width Seat	Depth (outside legs)	Depth Seat	Seat Height (highest)	Seat Height (lowest)	Product Weight	Carton Shipping Weight	Weight Limit
RTL12505	17"	21"	14"	2.25"	20.5"	16"	9 lbs.	10 lbs.	350 lbs.

bath safety

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM