



READY SET GO Deluxe Aluminum Shower Chair With Carry Bag

Item #	Product Description	Assembly	UOM
RTL12105KDR	Comes with Carry Bag	Knock Down	2/cs

- Easy, tool-free assembly of legs and back. (Figure A - Figure B)
- Drainage holes in seat and back reduce slipping.
- Aluminum frame is lightweight, durable and corrosion proof.
- New pinch-free cover allows for push pins to be depressed without pinching finger. The audible "click" will be heard when secure.
- Legs are height adjustable in 1" increments, and angled with suction style tips to provide additional stability.
- Support collar prevents leg movement and reduces rattle.
- Comes with a durable, nylon carry bag for easy and convenient transport. (Figure C)

specifications

Outside Legs	19" (W) x 18" (D)
Seat Dimensions	19.75" (W) x 12" (D)
Seat Height	14" - 18"
Weight	8 lbs.
Carton Shipping Weight	19 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

Shower Chair with Oval Back

Item #	UOM
RTL12605	1/ea

- Adds safety and support in the bath tub and shower.
- Easy, tool-free assembly and height adjustable legs.
- Corrosion-proof aluminum is both lightweight and durable.
- Non-slip feet.
- Built-in shower spray holder.
- Retail packaged.

specifications

Outside Legs	17.25" - 18"
Seat Dimensions	19.75" - 11.375"
Seat Height	14" - 21"
Weight	6.61 lbs.
Shipping Weight	8.38 lbs.
Weight Limit	300 lbs.
Warranty	Limited lifetime