

























Premium Plastic Raised, Regular/ Elongated Toilet Seat, with Lock

Item #	Product Description	UOM
12014	Without Arms	1/cs

- Designed for individuals who have difficulty sitting down or standing up from the toilet.
- · Wide opening in front and back for personal hygiene.
- · Lightweight and portable.
- Standard locking mechanism ensures safety. (Figure A)
- · Fits most elongated toilets.
- · No tools required for installation.
- Easy to clean.

Hinged Toilet Seat Riser

Item #	Product Description	UOM
RTL12607	Standard Seat	1/ea
RTL12608	Elongated Seat	1/ea

- Raises the height of the toilet 3.5" for individuals that require assistance sitting on and standing from the toilet.
- Hinge design allows riser to be moved out of the way when not in use. (Figure A)
- Heavy-duty molded plastic construction provides strength and durability.
- Lightweight, portable, and easily attaches to toilet bowl.
- Easy to install on most standard toilets.
- Two models to accommodate either standard or elongated (Figure B) toilet bowls.

specifications	
Seat Dimensions	17" (W) x 16.5" (D) x 4.5" (H)
Weight	4 lbs.
Carton Shipping Weight	6 lbs.
Weight Limit	300 lbs.
Warranty	Limited Lifetime

specifications	
Seat Dimensions (12607)	13.5" (W) x 17.5" (D) x 3.5" (H)
Seat Dimensions (12608)	14" (W) x 19.25 (D) x 3.5" (H)
Weight (12607 12608)	5.07 lbs. 5.62 lbs.
Weight Limit	250 lbs.
Warranty	Limited Lifetime