



PainAway Wireless TENS

Item

RTLAGF-900

UOM

1/ea

- Can be used for chronic and acute pain on most areas of the body.
- Wireless design with rechargeable remote and receivers – no batteries to purchase.
- Sleek LCD Remote controls up to 3 wireless receivers with 'snap connectors' to electrodes.
- Long lasting electrode pads for up to 30 uses.
- Dual Channel system delivers twice the power as conventional TENS units.
- Provides constant current throughout each treatment.
- Amplitude range: maximum: 50V-90V; minimum 0-10V 1KΩ load).
- Frequency: 1.2-60Hz $\pm 15\%$).
- Pulse: 100 μ S $\pm 30\%$).
- Comes complete with LCD Remote Control, Two Wireless Receivers, One Large Long Lasting Back Electrode pad, Two 2"x3" Long Lasting Electrode Pads, Recharging Kit, Quick Start Guide and Instruction Manual.
- Electrode refills available for separate purchase.
- Warranty: 2-Year, Remote & Receivers.