



Moist-Dry Heating Pad, Standard

Moist-Dry Heating Pad, Large

Item #	Product Size	UOM
RTLAF-HP-STD	Standard	6/cs

Item #	Product Size	UOM
RTLAF-HP-LG	Large	6/cs

- Delivers soothing moist heat to a variety of areas. (Figure A)
- Heat penetrates tissue quickly and safely.
- Cloth cover is machine washable.
- Convenient slide on/off switch with indicator light.
- Absorbent sponge included for moist heat application.
- Size: 11.5" x 13.5".
- Retail Packaged.

- Delivers soothing moist heat to a variety of areas. (Figure A)
- Heat penetrates tissue quickly and safely.
- Cloth cover is machine washable.
- 3 heat settings with slide switch and auto off provides comfort, safety and versatility.
- Self-adhesive belt ensures proper placement.
- Absorbent sponge included for moist heat application.
- Size: 12" x 24".
- Retail Packaged.

specifications

Warranty Limited 3 Years

specifications

Warranty Limited 5 Years