



Spitfire Scout 4 Compact Travel Scooter

- Quick and easy disassembly. (Figure A)
- Interchangeable color panels hide scuffs and scratches.
- Armrests are padded and adjustable.
- Quick connect batteries.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Delta tiller.
- Easy freewheel operation

| specifications | | | |
|----------------------|------------------------|--|-------------------|
| HCPCS | K0800 | Charger | 2A Offboard |
| Weight Capacity | 300 lbs. | Brakes | Electromagnetic |
| Top Speed | 4.25 mph | Front Wheels | 8" x 2" Flat-free |
| Maximum Range | 9 miles | Rear Wheels | 8" x 2" Flat-free |
| Turning Radius | 53.75" | Base Weight | 58 lbs. |
| Climbing Angle | 6° | Battery Weight (pair) | 20 lbs. |
| Ground Clearance | 2.5" | Seat Weight | 16 lbs. |
| Dimensions | 42.5" (L) x 21.25" (W) | Total Weight | 94 lbs. |
| Floor-to-Seat Height | 17" - 19" | Heaviest Piece Weight (Front Section) | 39 lbs. |
| Seat Dimensions | 16.5" (W) x 16" (D) | Warranty on Frame | Lifetime |
| Controller | CTE 45A | Warranty on Electronic Controller/Drive Train Components 14 months | |
| Motor | 24V x 270W x 4700 rpm | Warranty on Batteries | 6 months |
| Batteries (pair) | 12V x 12AH | | |