

ActiveCare
drive
 Mobility
SPITFIRE EX 1320

NEW!



**INTERCHANGEABLE
 COLOR PANELS
 (RED, BLUE & SILVER)
 INCLUDED!**



Midnight Blue Metallic Red Silver

**Spitfire EX 1320
 3 Wheel Compact Size Scooter**

- 1P** SPITFIRE132016FS12 Metallic Red, Folding Seat, 12Ah Battery, 1/cs
- 1P** SPITFIRE132016FS21 Midnight Blue, Folding Seat, 21Ah Battery, 1/cs

FEATURES:

- Scooter can be easily disassembled into 5 pieces for transport or storage.
- 9 mile maximum range with standard 12AH batteries.
- Solid, non-marking, flat-free tires for worry-free travel.
- Adjustable angle tiller allows the client to position the controller into the most ideal driving position.
- Seat features 180 degree swivel and fold-down backrest.
- Removable flip-up armrests have adjustable height, width and angle for a broad range of adjustments. These can easily be changed and adjusted to fit the needs of the client.
- Anti-tip wheels for greater safety.
- Seat back folds completely flat, which comes in handy when transporting or storing the scooter.
- Standard plastic front basket
- 300 lbs. Weight Capacity

**POV GROUP 1
 HCPCS
 CODE:
 K0800**

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

1P = 1-Pack Available

To order call toll free: 877.224.0946

POWERED MOBILITY

ActiveCare
 Mobility
drive