

## Features



- Disassembles quickly and easily, without tools, into 5 pieces, the heaviest weighing only 28 lbs.
- Highest-in-class 300 lb. weight capacity
- 9 mile maximum range (15 miles with optional 21AH batteries)
- Adjustable speed up to a maximum of 5 mph
- Adjustable angle tiller with fingertip throttle control
- Includes three sets of interchangeable body panels in red, blue and silver
- Quick connect battery pack
- Anti-tip wheels for greater safety
- Flip-up armrests with adjustable width and angle
- Anti-tip wheels for greater safety
- Comfortable, 17" wide swivel seat with fold-down, padded backrest
- Includes plastic front basket with removable liner bag
- Flat-free, non-marking tires
- Batteries included



## Included Color Panels

Metallic Red

Blue

Silver

## Accessories

- Oval Bag
- Carry All Bag
- Dust Cover

## Specifications

<b>Max Weight Capacity</b>	<b>300 lbs.</b>
<b>Max Speed</b>	5 mph
<b>Max Range</b>	9 miles (12AH)   15 miles (21AH)
<b>Max Grade/Incline</b>	12°
<b>Turning Radius</b>	48"
<b>Ground Clearance</b>	2"
<b>Dimensions</b>	
Overall Length	40"
Overall Width	20"
Floor to Seat Height	19"-22"
<b>Seating Dimensions</b>	
Seat Size	17"(w) x 16"(d)
Back Height	14"
Armrest to Seat Height	7"
<b>Power</b>	
Motor	24V x 270W x 4700rpm
Controller	Dynamic R50A
Batteries	Two DC12V x 12AH
Battery Charger	2 Amp, Offboard

<b>Brakes</b>	Electromagnetic
<b>Freewheel Mode</b>	Yes
<b>Tires (Flat-free)</b>	
Front	7" x 2"
Rear	8" x 2"
<b>Weights</b>	
Without Batteries	92 lbs.
With Batteries	111 lbs. (12AH) 121 lbs. (21AH)
Heaviest Piece Wt.	28 lbs. (rear section)

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

**ActiveCare**  
Medical

2 Harbison Way  
Columbia, SC 29212  
866.621.4933  
activecaremed.com

Information contained herein is correct at the time of publication. ActiveCare reserves the right to alter specifications without notice.