

## SILVER SPORT 1




### Silver Sport 1

Single Axle (K0001)  
(Has replaced Winnie and Saturn)

#### FEATURES:

- Powder coated silver vein steel frame.
- Padded armrests.
- Nylon upholstery is durable, lightweight, attractive and easy to clean.
- Urethane tires mounted on composite wheel provides durability and low maintenance.
- Comes with push to lock brakes.
- Carry pocket on backrest provides additional convenience.
- Comes with swing-away footrests or elevating legrests with tool free adjustable length front riggings.
- Plastic footplates.
- Color: Silver Vein.
- 3 Year Limited Warranty.



*Due to product improvements and changes, all images in this catalog are subject to change.*

*\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.*