

# BARIATRIC SENTRA EC HEAVY DUTY



**NOW  
AVAILABLE  
IN UP TO 24",  
AND 450 LB.  
WEIGHT  
CAPACITY.**

**BARIATRIC  
450 LB  
WEIGHT CAPACITY**

**NOW  
COMES WITH  
RIGGINGS IN  
CARTON!**

## **Bariatric Sentra EC Heavy Duty**

Deluxe Dual Axle (K0007)

### **FEATURES:**

- **Standard dual axle provides easy transition of seat to hemi-level.**
- **"Comes with Dual Cross Bar".**
- Carbon Steel frame with triple coated chrome for an attractive, chip-proof, maintainable finish.
- Double embossed vinyl upholstery is durable, light-weight, attractive and easy-to-clean.
- Urethane tires mounted on composite wheel provides durability and low maintenance.
- Aluminum footplates.
- **Comes with padded arm rests.**
- No maintenance, sealed, precision bearings.
- **8" front casters, adjustable in 3 positions.**
- Chrome handrims standard.
- Comes with push to lock wheel locks.
- Color: black upholstery
- Limited 3 Year Warranty.

*Due to product improvements and changes, all images in this catalog are subject to change.*

*\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.*