

BARIATRIC SENTRA EC HEAVY DUTY, EXTRA, EXTRA WIDE, (STANDARD)



Sentra EC Heavy Duty, Extra, Extra Wide (K0007)

BARIATRIC
700 LB
WEIGHT CAPACITY

FEATURES:

- Steel, spoke wheels are strong and maintenance free.
- 700 lb. Weight Capacity.
- Carbon steel frame with TRIPLE coated chrome for an attractive, chip-proof, maintainable finish.
- Reinforced steel gussets at all weight bearing points provide additional strength.
- Durable reinforced nylon upholstery.
- Extra large, padded armrests provide comfort and support.
- Precision sealed wheel bearings in front and rear ensure long lasting performance and reliability.
- Extra large heavy duty, aluminum footplates provide maximum support, are attractive and crack proof. (Must be ordered separately)
- Heavy duty, flat free 8" front casters are adjustable in 3 positions..
- Extra heavy duty front forks.
- Reinforced side frame, caster journals and dual cross bars.
- Comes with push to lock wheel locks.
- Color: Black upholstery, chrome frame.
- Limited Lifetime Warranty.

Due to product improvements and changes, all images in this catalog are subject to change.

**Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.*