



STDS1022

## Oblique Projection Handrim

STDS1022 1/ea

For Use with Sentra EC, Sentra Heavy Duty, Sentra Reclining, Viper Reclining, Viper, Viper Plus and Cirrus IV wheelchairs

### FEATURES:

- Provides user with additional "grip" when propelling a wheelchair.
- Easily attaches to wheel.
- Soft, rubber coating is comfortable and long lasting.



STDS1087

## Light Weight 8" Caster

STDS1087 1/ea

For use with Sentra Reclining, Viper Reclining, Viper, Viper Plus GT and Cruiser III wheelchairs

### FEATURES:

- Low profile, 8 spoke, nylon wheel
- 8" x 1" Pyramid Tire
- For use with all Drive wheelchairs.
- Limited 3 Year Warranty.



STDS840

## Flat-Free Wheel with Handrim

STDS840 1/bx

Includes Tire and Complete Wheel.

For use with Sentra Reclining, Viper Reclining, Viper, Viper GT and Cirrus IV wheelchairs

## Light Weight 24" Wheel

STDS1085 1/ea

For use with Sentra Reclining, Viper Reclining, Viper, Viper GT and Cirrus IV wheelchairs

### FEATURES:

- Custom polyurethane no-flat tire.
- Durable design provides shock absorption for greater comfort and reduces user fatigue.
- Composite hook-beak rim keeps tires on and prevents roll-offs.
- Integrated rim and hub allow years of trouble-free operation.
- 300 lb. Weight Limit.
- For use with all Drive wheelchairs.
- Limited 3 Year Warranty.



STDS1085

## One Arm Drive

STDS4055R Right, 1/ea

For use with 18" Viper and Viper Plus GT wheelchairs

STDS4055L Left, 1/ea

For use with 18" Viper and Viper Plus GT wheelchairs

STDS4058R Right, 1/ea

For use with 18" Sentra EC wheelchair

STDS4058L Left, 1/ea

For use with 18" Sentra EC wheelchair



STDS4055

### FEATURES:

- Allows wheelchair to be propelled with one arm.

*Due to product improvements and changes, all images in this catalog are subject to change.*

*\*Important Notice: Item numbers that display an asterisk (\*) reflect an item that will be discontinued or package change.*

wheelchairs

To order call toll free: 877.224.0946

drive

WWW.DRIVEMEDICAL.COM