

**NEW & IMPROVED
WITH STRONGER
HINGED ARM!**



STDS1040S



**INSTALLATION
AND REMOVAL
IS EASY!**



Universal Cup Holder

STDS1040S (Has replaced STDS1040), 1/ea

FEATURES:

- Easily attaches to wheelchairs, rollators or walkers with tubing diameter from 5/8" to 1".
- Special clamp insert can be added or removed as necessary to accommodate different size tubing.
- Made with durable, lightweight plastic.
- Installation and removal is easy with adjustment knob and clamp.
- Dimensions: 3" (W) x 3" (H).
- Limited 1 year warranty.



STDS5050

Wheelchair Tray

STDS5050 Easy To Clean Surface,
Tray Size: 24" x 20" x 1/2", 1/cs
Accommodates Wheelchairs
16" - 20" Seat Width



Grommets allow for use with 16", 18"
and 20" wheelchairs.



STDS835

Chart Carry Pocket

STDS835 Nylon, 1/ea
For use with 16", 18" and 20"
Wheelchairs with Nylon Upholstery

STDS836 Vinyl, 1/ea
For Use with 16", 18" and 20"
Wheelchairs with Vinyl Upholstery

*Due to product improvements and changes,
all images in this catalog are subject to change.*

**Important Notice: Item numbers that display an asterisk (*) reflect an item
that will be discontinued or package change.*

wheelchairs

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM