



STDS810

Seat Rail Extension

Adds 2" to depth of wheelchair seat.

- STDS810 Black, 1pr/bx (shown)
For use with Viper and Viper GT
(16" and 18" Seat Widths Only)
- STDS811 Chrome, 1pr/bx
For use with Sentra EC 20"
- STDS813 Black, 1pr/bx
For use with Cirrus IV
(16" and 18" seat widths only)



STDS860

Universal Forward Stabilizers

- STDS860 Chrome, 1pr/bx
Universal spacing sleeve can be added or removed to accommodate a secure fit with Sentra EC 20". Sentra HD 20", Sentra Reclining, Cruiser III, Cirrus IV, Viper, Viper Plus GT, Viper Plus Reclining and most other leading manufacturers wheelchairs
- STDS862 Chrome, 1pr/bx
For use with 20", 22" and 24" Sentra EC wheelchairs



STDS850

Seat Belt

- STDS850 Auto Type, Black, 1/bx (shown)
- STDS851 Velcro Type Closure, Black, 1/bx
For use with all Drive and other leading manufacturers wheelchairs.
Can accommodate waist sizes up to 48 inches.



STDS851

Bariatric Seat Belts

- STDS855 Auto Type, Black, 60", 1/bx
- STDS856 Velcro Type, Black, 60", 1/bx
For use with all Drive and other leading manufacturers wheelchairs.
Can accommodate waist sizes up to 60 inches.



STDS855

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

wheelchairs

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM