

TRAVELITE TRANSPORT CHAIR



COMES COMPLETE WITH READY-SET-GO CARRY BAG!



LIGHT-WEIGHT, EASY TO CLEAN UPHOLSTERY!



TRANSPORT CHAIR FOLDS INTO COMPACT UNIT FOR STORAGE AND TRANSPORTATION!



Travelite Transport Chair

TC005 Chair in a Bag, 1/cs

FEATURES:

- The Travelite Chair in a Bag is a lightweight yet sturdy aluminum transport chair in gunmetal grey.
- It is easily transportable and comes with a compact carry bag making it ideal for traveling.
- It is ideal for use indoors and outdoors.
- Folding backrest.
- 18" square nylon seat canvas with lift up straps.
- Comes with parking brakes
- Comes with padded, flip up desk arms and height adjustable footrests.
- 8" solid tires.
- Weight capacity 250 lbs.
- Limited Lifetime Warranty

SPECIFICATIONS

	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length (with riggings)	Product Weight	Carton Shipping Weight	Weight Limit
18"	22.5"	11.5"	37.5"	19"	17.5"	18"	18"	33"	20.5 lbs.	23.2 lbs.	250 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

transport chairs

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM