



TR39E-SV

WEIGHS ONLY 26 lbs.!



CONVENIENT CUP HOLDER
CAN EASILY ATTACH TO CHAIR!
SOLD SEPARATELY MODEL # STDS1040S

**NEW SILVER
VEIN FINISH!**



**CARRY BAG
AVAILABLE**

SEE PAGE WC 278

Steel Transport Chair

TR37E-SV	17" Silver Vein Frame and Black Upholstery, (Has replaced TR27E-SV), 1/cs
TR39E-SV	19" Silver Vein Frame and Black Upholstery, (Has replaced TR29E-SV), 1/cs
TR17SV-DDA	17" Silver Vein Frame, Black Upholstery w/Removable Desk/Arm, (Has replaced TR17-DDA), 1/cs
TR19SV-DDA	19" Silver Vein Frame, Black Upholstery w/Removable Desk/Arm, (Has replaced TR19-DDA), 1/cs
LELR-TF	Swing-away Elevated Legrest, Tool Free, (Optional front riggings), 1pr/cs

FEATURES:

- Attractive blue or silver vein finish.
- Weighs 26 lbs.
- Available in 17" or 19" seat widths.
- Back folds down for storage and transport.
- Swing-away footrests come standard.
- **Comes with seat belt for added safety.**
- Padded armrests provide additional comfort.
- 8" casters in front and rear.
- Comes with attractive nylon upholstery that is durable and easy to clean.
- Standard rear wheel locks.
- Comes with upholstery handles to fold chair.
- Limited Lifetime Warranty.



**BACK FOLDS DOWN FOR
EASY FOLDING AND STORAGE.**

SPECIFICATIONS

Item #	Overall Width (open)	Overall Width (closed)	Overall Height	Seat To-Floor Height	Back Height	Width Seat	Depth Seat	Overall Length w/Rigging	Product Weight	Carton Shipping Weight	Weight Limit
TR37E-SV	18"	8"	39"	19"	18"	17"	16"	33"	26 lbs.	31 lbs.	250 lbs.
TR17SV-DDA	20.5"	10.75"	37.5"	19"	17"	17"	16"	33"	29 lbs.	34 lbs.	300 lbs.
TR39E-SV	20"	8"	39"	19"	18"	19"	16"	33"	26 lbs.	31 lbs.	250 lbs.
TR19SV-DDA	23"	10.75"	37.5"	19"	17"	19"	16"	33"	29 lbs.	34 lbs.	300 lbs.

Due to product improvements and changes, all images in this catalog are subject to change.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.