



COMPARE & SAVE!

A less expensive alternative
to the Convaid® Cruiser®!

Trotter Mobility Chair

NEW!

The Trotter is a lightweight portable positioning chair that provides transport solutions for parents and children. A wide range of options provide custom positioning for your individual needs and its simple adjustments allow the Trotter to grow with your child.

TR 1200	Trotter, 12", 1/cs
TR 1400	Trotter, 14", 1/cs
TR 1600	Trotter, 16", 1/cs
TR 1800	Trotter, 18", 1/cs

FEATURES:

- Lightweight folding frame with closure strap that holds the chair shut when folded.
- Black removable and washable fabric.
- Crash tested for transit. The Trotter has passed a dynamic crash test done with a 130 lb. front-facing occupant tested at 30 mph.
- Adjustable seat angle - 15° and 22°.
- Adjustable back angle - 85°, 90° and 95°.
- 8" front casters and 10" rear wheels with toggle lock.
- Adjustable seat depth.
- Curb assist lever facilitates curbside navigation.
- Positioning pelvic belt and five point "H" harness with padded strap covers. The reversible padded strap covers increase comfort and one side of cover is imitation lamb's wool and the other is velour.
- Swing-away footrest are independently adjustable and removable.
- Two Year Limited Warranty.



**CURB
ASSIST
LEVER**



ITEM #	SEAT WIDTH	SEAT DEPTH	BACK HEIGHT	SEAT TO FOOTPLATE	SEAT TO FLOOR	PRODUCT WEIGHT WITHOUT FRONT RIGGING	CARTON SHIPPING WEIGHT	WEIGHT LIMIT STANDARD	WEIGHT LIMIT WITH TRANSIT OPTION
TR 1200	11.75"	11.5" - 13.5"	20"	8" - 22"	22"	19 lbs.	28 lbs.	75 lbs.	70 lbs.
TR 1400	13.75"	13" - 15"	23"	9" - 22"	22"	21 lbs.	30 lbs.	100 lbs.	100 lbs.
TR 1600	15.5"	13.5" - 15.5"	24"	10" - 22.5"	22"	24.5 lbs.	33.5 lbs.	170 lbs.	170 lbs.
TR 1800	17.5"	15" - 17"	25"	11" - 23"	22.5"	30 lbs.	39 lbs.	250 lbs.	170 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.



ACCESSORIES:

All accessories are universal and fit all four Trotter sizes.

Headrest Extension

TR 8021 Headrest Extension, 1/bx

Black nylon fabric Headrest Extension is width and height adjustable. It provides additional head support for taller users. Used on bus transit models to help prevent head whiplash during rear end collision.

Padded Headrest Wing

TR 8022 Padded Headrest Wing, 1/bx

Height adjustable headrest wing provides lateral support for the head.

Utility Bag

TR 8023 Utility Bag, 1/bx


Zippered bag hooks on to the back of the Trotter mobility chair.

Full Torso Vest

TR 8025 Full Torso Vest, 1/bx

Adjustable vest provides support and maintains proper body posture and alignment. It replaces the "H" Harness.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.



mobility

mobility chair

trotter

Wenzel *ite Re/hab*

To order call toll free: 877.224.0946

drive

ACCESSORIES:

Canopy

TR 8026 Canopy, 1/bx

Height adjustable and removable canopy provides protection from sun or rain and includes a clear flap-covered observation panel.

Upper Extremity Support Tray

TR 8024 Clear Tray, 1/bx

Trays feature a safety lip and are easy to clean. Provides trunk stabilization and are great for feeding.

Lateral Support & Scoli Strap

TR 8027 Lateral Support and Scoli Strap, 1/bx

Adjustable and removable strap provides lateral support. It can be pulled to one side for Scoliosis correction.

Foot & Ankle Positioner

TR 8028 Foot & Ankle Positioner, 1pr/bx

Foot and ankle straps provide adjustable support and positioning.

Bus Transit Tie-Downs

TR 8029 Bus Transit Option, 1 Set of 4/bx

Four tie-down brackets allow the Trotter to be used for transit. (Must also order headrest extension & foot/ankle positioner to be in compliance with WC-19 Safety Standard Testing)

Seat Rail Extension Kit

TR 8031 Seat Rail Extension Kit, 1pr/bx

The seat rail extension kit adds an additional 3" of seat depth.

Colored Upholstery

- TR 12SB-B Blue Upholstery for TR 1200, 1/bx
- TR 12SB-R Rose Upholstery for TR 1200, 1/bx
- TR 14SB-B Blue Upholstery for TR 1400, 1/bx
- TR 14SB-R Rose Upholstery for TR 1400, 1/bx
- TR 16SB-B Blue Upholstery for TR 1600, 1/bx
- TR 16SB-R Rose Upholstery for TR 1600, 1/bx
- TR 18SB-B Blue Upholstery for TR 1800, 1/bx
- TR 18SB-R Rose Upholstery for TR 1800, 1/bx

Choice of colored upholstery easily replaces standard black upholstery.



Foot & Ankle Positioner



Scoli Strap



Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.