



Item #	Seat	UOM
VENTURA318CS	18" Captain's Seat	1/cs
VENTURA320CS	20" Captain's Seat	1/cs



Ventura 3 DLX

Midsized Scooter

- Includes interchangeable color panels in Red and Blue. (Figure A)
- 36 AH batteries for extended range.
- Armrests are padded and adjustable.
- Floorboard with 23" of foot space.
- Ergonomic throttle control.
- Height-adjustable swivel seat features fold-down backrest.
- Seat position can be adjusted forward or backward.
- Rear suspension.
- Programmable controller.
- Seat back reclines and features adjustable headrest.
- Available with 18" or 20" width Captain's Seat.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Headlight.
- Rearview mirrors.
- Turn signals/hazard lights.
- Delta tiller.
- Easy freewheel operation.

specifications

HPCPS	N/A	Batteries (pair)	12V x 36AH
Weight Capacity	350 lbs.	Charger	4A Offboard
Top Speed	5 mph	Brakes	Electromagnetic
Maximum Range	22 miles	Front Wheels	10" x 3" Flat-free
Turning Radius	53"	Rear Wheels	10" x 3" Flat-free
Climbing Angle	7.5°	Base Weight	88 lbs.
Ground Clearance	2.5"	Battery Weight (pair)	46 lbs.
Dimensions	48" (L) x 25" (W)	Seat Weight	40 lbs.
Floor-to-Seat Height	20.5" - 23.5"	Total Weight	174 lbs.
Seat Type	Captain's	Warranty on Frame	Lifetime
Seat Dimensions	18" (W) x 18" (D) 20" (W) x 18" (D)	Warranty on Electronic Controller/Drive Train Components	14 Months
Motor	24V x 350W x 5100 rpm	Warranty on Batteries	6 Months