



Item #	Seat	UOM
VENTURA318FS	18" Folding Seat	1/cs



HPCPS  
**K0807**



## Ventura 3

### Midsize Scooter

- Includes interchangeable color panels in Red and Blue. (Figure A)
- Armrests are padded and adjustable.
- Floorboard with 23" of foot space.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.
- Flat-free, non-marking tires.
- Seat position can be adjusted forward or backward.
- Rear suspension.
- Programmable controller.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Headlight.
- Delta tiller.
- Easy freewheel operation.

#### specifications

HPCPS	K0807	Batteries (pair)	12V x 20AH
Weight Capacity	350 lbs.	Charger	4A Offboard
Top Speed	5 mph	Brakes	Electromagnetic
Maximum Range	15 miles	Front Wheels	10" x 3" Flat-free
Turning Radius	53"	Rear Wheels	10" x 3" Flat-free
Climbing Angle	7.5°	Base Weight	88 lbs.
Ground Clearance	2.5"	Battery Weight (pair)	28 lbs.
Dimensions	48" (L) x 25" (W)	Seat Weight	26 lbs.
Floor-to-Seat Height	21.5" - 23.5"	Total Weight	142 lbs.
Seat Type	Folding	Warranty on Frame	Lifetime
Seat Dimensions	18" (W) x 17.5" (D)	Warranty on Electronic Controller/Drive Train Components	14 Months
Motor	24V x 350W x 5100 rpm	Warranty on Batteries	6 Months