



| Item # | Seat | UOM |
|--------------|--------------------|------|
| VENTURA418CS | 18" Captain's Seat | 1/cs |
| VENTURA420CS | 20" Captain's Seat | 1/cs |



Ventura 4 DLX

Midsize Scooter

- Includes interchangeable color panels in Red and Blue. (Figure A)
- Armrests are padded and adjustable.
- 36AH batteries for extended range.
- Ergonomic throttle control.
- Height-adjustable swivel seat features fold-down backrest.
- Seat back reclines and features adjustable headrest.
- Seat position can be adjusted forward and backward.
- Front and rear suspension.
- Programmable controller.
- Available with 18" or 20" width Captain's Seat.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Headlight.
- Rearview mirrors.
- Turn signals/hazard lights.
- Delta tiller.
- Easy freewheel operation.

specifications

| | | | |
|----------------------|---------------------------------------|--|--------------------|
| HCPCS | N/A | Batteries (pair) | 12V x 36AH |
| Weight Capacity | 350 lbs. | Charger | 4A Offboard |
| Top Speed | 5 mph | Brakes | Electromagnetic |
| Maximum Range | 22 miles | Front Wheels | 10" x 3" Flat-free |
| Turning Radius | 54" | Rear Wheels | 10" x 3" Flat-free |
| Climbing Angle | 7.5° | Base Weight | 99 lbs. |
| Ground Clearance | 2.5" | Battery Weight (pair) | 46 lbs. |
| Dimensions | 48" (L) x 25" (W) | Seat Weight | 40 lbs. |
| Floor-to-Seat Height | 22" - 24" | Total Weight | 185 lbs. |
| Seat Dimensions | 18" (W) x 18" (D) 20" (W) x 18" (D) | Warranty on Frame | Lifetime |
| Seat Type | Captain's | Warranty on Electronic Controller/Drive Train Components | 14 Months |
| Motor | 24V x 350W x 5100 rpm | Warranty on Batteries | 6 Months |