



Item #	Seat	UOM
VENTURA418FS	18" Folding Seat	1/cs



HCPCS
K0807



Ventura 4

Midsize Scooter

- Includes interchangeable color panels in Red and Blue. (Figure A)
- Armrests are padded and adjustable.
- Ergonomic throttle control.
- Height-adjustable swivel seat with fold-down backrest.
- Flat-free, non-marking tires.
- Large, plastic carry basket.
- Seat position can be adjusted forward or backward.
- Front and rear suspension.
- Programmable controller.
- Easy-to-adjust tiller.
- Anti-tip wheels.
- Headlight.
- Delta tiller.
- Easy freewheel operation.

specifications

HCPCS	K0807	Batteries (pair)	12V x 20AH
Weight Capacity	350 lbs.	Charger	4A Offboard
Top Speed	4 mph	Brakes	Electromagnetic
Maximum Range	15 miles	Rear Wheels	10" x 3" Flat-free
Turning Radius	54"	Base Weight	99 lbs.
Climbing Angle	7.5°	Battery Weight (pair)	28 lbs.
Ground Clearance	2.5"	Seat Weight	26 lbs.
Dimensions	48" (L) x 25" (W)	Total Weight	153 lbs.
Floor-to-Seat Height	22" - 24"	Warranty on Frame	Lifetime
Seat Type	Folding	Warranty on Electronic Controller/Drive Train Components	14 Months
Seat Dimensions	18" (W) x 17.5" (D)	Warranty on Batteries	6 Months
Motor	24V x 350W x 5100 rpm		