



Wallaby Pediatric Wheelchair

Item #	Product Seat Dimensions	Frame Color	Product Overall Width Open	Product Overall Width Closed	Product Overall Length (with riggings)	Product Length	Product HCPCS	Product Weight Capacity	UOM
WB1200-2GFR	12" (W) x 12" (D)	Red	20"	12"	37"	28.5"	E1238	150 lbs.	1/ea
WB1400-2GJB	14" (W) x 14" (D)	Blue	22"	12"	39"	30.5"	E1238	150 lbs.	1/ea

- Standard Transit Tie Down Brackets.
- WC19 Compliant when used with WB 8021 for users between 50 and 80 lbs.
- Height adjustable push handles.
- Flip back padded desk arms.
- Pelvic belt standard.
- Push to lock wheel locks.
- 22" Rear wheels with 6" front casters.
- Swingaway footrests with tool free height adjustment and calf strap.

specifications

Seat to Floor Height	19"
Backrest Height	14"
Armrest Height	7"
Weight Capacity	150 lbs.
Weight Capacity when used in transport	50-80 lbs.
Product Weight with footrests	37 lbs.

accessories

WB 8021	Headrest Extension	1/bx
WB ELR	Elevating Legrests for Wallaby	1/pr
STDS833N	Anti-Tippers With Wheels Black	1/pr