



Features

- Lightweight steel frame can be **folded** to a compact 13" width for easy portability and storage.
- General foam **seat cushion** included.
- **Oversized padded, flip-up armrests** allow for easy transfers and access to user. Armrests are also height adjustable.
- **Adjustable length controller mount** can be affixed to either the left or right armrest and can be fitted with an optional **swingaway arm** for additional flexibility.
- Standard **swingaway foot rests** can be flipped up out of the user's way and height adjusted.
- Soft plastic **skirt guards** eliminate sharp edges against client and clothing.
- Nylon seat and back **upholstery** offers more breathable seating surface, reducing heat and moisture buildup.
- Standard **calf strap** adds additional support for legs.
- Fold-down **backrest** makes it easier for storage and transport, plus it offers additional access to user for an attendant or caregiver.
- Rear **anti-tip wheels** for added safety.
- Seat back **storage pouch**, handy storage for any take-along items for user or caregiver.
- Chair can be fitted with optional **Elevating Leg Rests** allowing the user to travel comfortably with their legs raised—a necessity for those with some medical conditions, such as edema.
- Standard **positioning belt** for added user security.
- Manual **wheel locks** offer additional security and safety.
- **Programmable controller** with common connectors is easy to service and setup and can be personalized for individuals with special requirements.
- **Flat-free, non-marking tires** offer worry-free travel.
- Motor assemblies offer easy **freewheel operation**, allowing the chair to be pushed, which can be helpful for transport.

WildcatTM

Folding Power Wheelchair

ActiveCare
Medical

Wildcat™ Folding Power Wheelchair

Specifications

HCPCS K0822

Weight Capacity 300 lbs.

Maximum Speed 4 mph

Maximum Range 20 miles

Turning Radius 32"

Ground Clearance 4"

Overall Length 43"

Overall Width 26" (18" Seat)
28" (20" Seat)

Floor to Seat Height 19"

Controller Dynamic A50A

Motors Two 24V x 200W x 3900rpm

Batteries Two 12V x 35AH (UI)

Battery Charger 5 Amp, Offboard

Brakes Electromagnetic

Caster Wheels 8" x 2" Solid

Drive Wheels 12.5" x 2.25" Flat-Free

Base Weight 92 lbs.

Battery Weight (pair) 48 lbs.

Speed and range will vary according to user weight, terrain, battery charge/condition and tire pressure.

Foldable for Easy Transport



The Wildcat disassembles and folds to a narrow 13" width for storage or transport. Simply disconnect and remove the two battery cases, remove the seat cushion and pull the chair's sides toward each other. You can also fold the backrest down to make it even more compact.

Options

Seat Options

- ☐ 18" (w) x 16" (d) x 19" (h) ☐ 20" (w) x 16" (d) x 19" (h)

Front Riggings

- ☐ Swingaway Foot Rests Included
☐ Elevating Leg Rests Sold Separately

Controller Arm

- ☐ Adjustable Length Included
☐ Swingaway Sold Separately

Armrests

- ☐ Flip-Up, Full Length, w/ Adj. Height & Width Included

Accessories

- ☐ Armrest Bag

Color Options



Black



Midnight Blue



Neon Yellow



Rose Petal Pink

ActiveCare
Medical

2 Harbison Way | Columbia, SC 29212
866.621.4933 | activecaremedical.com