



Seat2GO

Positioning Seat

PORTABLE • COMFORTABLE • CONVENIENT

Take it with you everywhere you go!!

A lightweight soft padded positioning seat for children that need support while sitting on the floor, at the table or at play. It is simple to use in a home, school or when traveling

CRS 2000 Small, 1/cs
CRS 3000 Medium, 1/cs

FEATURES:

- Contoured seat prevents child from sliding forward and improves seating posture.
- 95° seat to back angle.
- Double strap trunk support and belt.
- Two sets of straps secure Seat2Go to any chair. When not in use, the seat straps can be stored under the seat behind the clear flap and the back straps in the back pocket.
- Can be used as an insert for a single or double stroller.
- 1 Year Limited Warranty.



WEIGHS ONLY 6 LBS!



ITEM #	WIDTH	DEPTH	HEIGHT OF BACK	WEIGHT LIMIT	PRODUCT WEIGHT	SHIP WEIGHT
CRS 2000	8"	9"	12"	50 lbs.	5 lbs.	8 lbs.
CRS 3000	12"	12"	14"	75 lbs.	7 lbs.	11 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.

positioning

seating

Wenzel[®] Lite Re/hab

To order call toll free: 877.224.0946

drive
WWW.DRIVEMEDICAL.COM