

NEW!

**AUTO-ADJUSTING
CONSERVER SENSES
BREATHING PATTERN!**



The SmartDose® Mini Oxygen Conserving will sense patient breath rate and turn oxygen flow up and down as needed, within seconds.

Auto-adjusting conservers are the future of oxygen therapy, optimizing both patient saturations and maximizing conservation. Satisfy your patients' and referrals' wishes while decreasing your delivery costs, all with one unit.



SmartDose® Mini, Electronic Oxygen Conserving

1P CTOX-MNO₂ 1/ea

FEATURES:

- The SmartDose® Mini Oxygen conserving will sense patient breathing pattern and auto-adjust oxygen flow rate up or down as needed, within seconds
- Uniform pulse dose increments 1, 2, 3, 4, 5 plus CF preset @ 2 LPM
- Achieves maximum conservation efficiency by automatically adjusting the dose down at lower activity levels
- Operating pressure range: 500 - 3000 PSI
- Weighs just 15.5 ounces with batteries
- Delivers a quiet, gentle dose of oxygen using an efficient, compact design
- Can use 2 AA Alkaline batteries for a duration of about one year for most patients
- 2 Year Limited Warranty

ACCESSORIES:

OP-150-800	Oxygen Cylinder Shoulder Carry Bag, Fits M4, M6, M7, and M9, 1/ea
OP-150T	Oxygen Cylinder Horizontal Carry Bag, Fits M4, M6, M7, and M9, 1/ea
18102	Oxygen Cylinder Shoulder Carry Bag, Fits up to D Cylinders, 1/ea
Soft 207	Non-Kinking, 7' Cannula, 50/cs
Soft 007	Curved, 7' Cannula, 50/cs
Seal100	Replacement O-Ring Washer Seal, 10/bag

1P = 1-Pack Available

respiratory

To order call toll free: 877.224.0946

drive